

## CONNECTION CABLES



CE

Cod. **MNC12PDINS**  
CONNECTION CABLE 120 cm long with T.P. connector one side and a 1mm socket on the other side . One set 10 pcs.

Cod. **MNC12DDINS**  
CONNECTION CABLE 120 cm long with T.P. connector one side and a 1,5mm Pin on the other side. One set 10 pcs

Cod. **MNC12MEDDINS**  
CONNECTION CABLE 120 cm long with T.P. connector one side and a 0,7mm socket on the other side. One set 10 pcs




Cod. **JUMP2626**  
Jumper lead adaptor

CE


Cod. **1003SNAP71**  
Flexible shielded cable with 3 SNAP connector one side and a 5 Din connector the other side

Cod. **1003SNAP26**  
Flexible shielded cable with 3 SNAP connector one side and a 3 T.P. sockets the other side



Cod. **1002MED71**  
Flexible shielded cable with 2 0,7 mm socket an other side 5pole din

Cod. **1002MED26**  
Flexible shielded cable with 2 0,7 mm socket an other side 3 TP connectors




Cod. **1003DIN71**  
Flexible shielded cable with 3 male T.P. one side and a 5 Din connector the other side

Cod. **1003DIN26**  
Flexible shielded cable with flexible 3 male T.P. one side and a 3 T.P. sockets the other side




Cod. **1003P1N71**  
flexible shielded cable with three 1 mm socket one side and 5 Din pole connector the other side

Cod. **1003P1N26**  
flexible shielded cable with three 1 mm socket one side and 3 T.P. connectors the other side



Cod. **1003CROC71**  
cable with 3 crocodile connector one side an other side 5pole din

Cod. **1003CROC26**  
cable with 3 crocodile connector one side an other side 3 TP connectors



CE

## Silver/silver chloride disposable electrodes with adhesive and

### conductive gel for neurological application

The wire is melted to the electrode; this prevents the cable from becoming detached. The entire surface of the AgAgCL is conductive. The solid gel keeps the electrode well in place. The coloured lead wires are available with T.P sockets, snap connectors, tab connectors, 1 mm pin connectors or 0.7 mm pin connectors.


#### 15x20mm disposable electrode. 6 pcs box.



CE

Cod. DENIS01526 with 15cm long lead wire and Touch proof connector  
 Cod. DENIS10026 with 100cm long lead wire and Touch proof connector  
 Cod. DENIS15026 with 150cm long lead wire and Touch proof connector

#### 25x20mm disposable electrode. 6 pcs box.



CE

Cod. DENIL01526 with 15cm long lead wire and Touch proof connector  
 Cod. DENIL10026 with 100cm long lead wire and Touch proof connector  
 Cod. DENIL15026 with 150cm long lead wire and Touch proof connector

#### Disposable electrode with 8cm long lead wire and 0,7mm pin connector. 12 pcs box.



CE

Cod. DENIS00807  
 15x20mm disposable electrode with 8cm long lead wire and 0,7mm pin connector

Cod. DENIL00807  
 25x20mm disposable electrode with 8cm long lead wire and 0,7mm pin connector

#### Disposable electrode with 8cm long lead wire and 1mm pin connector. 10 pcs box.




CE

Cod. DENIS00830  
 15x20mm disposable electrode with 8cm long lead wire and 1mm pin connector

Cod. DENIL00830  
 25x20mm disposable electrode with 8cm long lead wire and 1mm pin connector

#### Disposable electrode with STUD connector



CE

#### Disposable TAB electrode




CE


Cod. TABSM00000  
 Disposable TAB electrode

## Disposable and reusable SUBDERMAL needle electrodes

As an alternative to disc electrodes the needle electrodes are available straight or hook-shaped. Different lengths are available on request the sharpness is 15° specially made for a non traumatic insertion. A special care is dedicated during the production of the needle: the needle and the cable are moulded together to offer a comfortable and safety handgrip.



### DISPOSABLE ETO sterilized SUBDERMAL needle electrode . 25pcs box.



Cod. MN4013D10SM diam.0,40x13mm needle with 100cm long 5 colors set cable and safety connector  
Cod. MN4013D15SM diam.0,40x13mm needle with 150cm long 5 colors set cable and safety connector  
Cod. MN4013D25SM diam.0,40x13mm needle with 250cm long 5 colors set cable and safety connector



0476

### DISPOSABLE SUBDERMAL needle electrode with 8cm long cable and 1mm pin connector



0476

Cod. MN3512P150  
disposable ETO sterilized diam.0,35x12mm needles with 8cm long cable and 1mm pin connector

### DISPOSABLE SUBDERMAL hook needle electrode with 8cm long cable and 1mm pin connector



0476

Cod. MN3512H150  
disposable ETO sterilized diam.0,35x12mm hook needles with 8cm long cable and 1mm pin connector

### DISPOSABLE SUBDERMAL needle electrode with 1mm pin connector



0476

Cod. MN3512F110  
disposable ETO sterilized diam.0,35x12mm needle with 1mm pin connector