

Bipolar surface Stimulation Electrodes and Recording Bar electrodes



CE

This Stimulator is for surface stimulation with felts or steel-balls. There is a 30mm distance between the two stimulation points and can be either hand held or fixed with Velcro bands. The stimulator electrode with hand grip is sold in a kit which contains:

- Round felter tips
- Two stainless- steel stimulating electrodes
- A screwdriver to screw the metal electrodes in the stimulator
- Two Velcro straps 40 and 70 cm long

Cod. **STIM1571**
150cm long lead wire and 5-pole DIN45322 connector

Cod. **STIM1526**
150cm long lead wire and Touch proof connector

Cod. **SFERSTIM**
stain&steal stimulating electrode to be used with cod. STIMXXXX

Cod. **FELT0520**
disposable felter to be used cod. STIMXXXX 1 pack 50 pz



CE

Small bipolar shielded bar electrodes.

Recording and stimulating bar electrodes with two 6,5mm x 20mm felts. 25mm distance between the electrodes. The anode is marked to distinguish between plus and minus as well as different and indifferent. All recording electrodes are sold in a kit which contain the stimulator / recorder with Two Velcro strap cm40 and cm 70 long.

Cod. **STRCP01571**
150cm long lead wire and 5-pole DIN45322 connector

Cod. **STRCP01526**
150cm long lead wire and Touch proof connector

Cod. **FELT0720**
disposable felter to be used (cod.STRCPXXXX) 1 pack 50 pz



CE

Large bipolar shielded bar electrodes.

Recording and stimulating bar electrodes with two 6,5mm x 20mm felts. 40mm distance between the electrodes. The anode is marked to distinguish between plus and minus as well as different and indifferent. All recording electrodes are sold in a kit which contain the stimulator / recorder with Two Velcro strap cm40 and cm 70 long.


Cod. **STRCG01571**
150cm long lead wire and 5-pole DIN45322 connector

Cod. **STRCG01526**
150cm long lead wire and Touch proof connector

Cod. **FELT0720**
disposable felter to be used (cod.STRCGXXXX) 1 pack 50 pz



Nerve surface stimulator



Hand grip bar with 10mm diam. stainless steel disc electrodes. 30mm distance between the electrodes.



Cod. **BAIM0026** 100cm long lead wire and Touch proof connector
Cod. **BAIM0071** 100cm long lead wire and 5-pole DIN45322 connector



Bar with 10mm diam. stainless steel disc electrodes. 30mm distance between the electrodes.



Cod. **BARR0026** 100cm long lead wire and Touch proof connector
Cod. **BARR0071** 100cm long lead wire and 5-pole DIN45322 connector

Stimulator with metal electrodes, with interchangeable distance between electrodes (from 8 to 15-23-30mm)



Cod. **STIM0171**
100cm long lead wire and 5-pole DIN45322 connector

Cod. **STIM0126**
100cm long lead wire and Touch proof connector



A pair of conductive soft stainless-steel tissue rings electrodes



A pair of conductive, soft, stainless steel tissue ring electrodes made and a clever fastener which don't permit the electrodes to move during the stimulation.(patent pending)

Cod. **SRCOLP0126** A pair of conductive rings with 100cm long cable and safety connector Din 42802


Cod. **SRCOLP0171** A pair of conductive rings with 100cm long cable and 5-pole DIN 45322 connector

Cod. **SRCOLPS126** A pair of small conductive rings with 100cm long cable and safety connector Din 42802

Cod. **SRCOLPS171** A pair of small conductive rings with 100cm long cable and 5-pole DIN 45322 connector



A pair of thin STAINLEES STEEL RINGS electrodes




Cod. **SRRING0126**
a pair of stainless steel rings with 100cm long cable and Touch proof connector

Cod. **SRRING0171**
a pair of stainless steel rings with 100cm long cable and 5-pole DIN45322 connector

Cod. **SRRINGS126**
a pair of small stainless steel rings with 100cm long cable and Touch proof connector

Cod. **SRRINGS171**
a pair of small stainless steel rings with 100cm long cable and 5-pole DIN45322 connector





CE

Cod. **STIM/REC00** Set of reusable stimulating and recording electrodes. The set contains:

- 1 Cod. **100000NH** Reusable EMG connection shielded cable with 5-pole connector for RECORDING
- 2 Cod. **1002SCNHR** Reusable EMG connection cable with 2 Touch Proof connector for STIMULATION

- A Cod. **SRSOC0000** A pair of 1mm sockets moulded in a bipolar socket to fit the EMG needle cable
- B Cod. **SRDISC0000** A pair of AgAgCl disc electrodes moulded in a bipolar socket to fit the EMG needle cable
- C Cod. **SRCOLP0000** A pair of conductive soft stainless- steel tissue rings moulded in a bipolar socket to fit the EMG needle cable
- D Cod. **SRCLIP0000** A pair of insulated alligator clips moulded in a bipolar socket to fit the EMG needle cable
- E Cod. **SRRING0000** A pair of stainless-steel rings moulded in a bipolar socket to fit the EMG needle cable
- F Cod. **SRFLAT0000** A pair of AgAgCl flat electrodes moulded in a bipolar socket to fit the EMG needle cable
- G Cod. **SRSTUD0000** A pair of stud adapters moulded in a bipolar socket to fit the EMG needle cable

GROUND ELECTRODES

"PLUS" Ground electrodes

100cm long lead wire



This special ground electrode is made of a very flexible mesh of AISI 316L stainless-steel with 0 Ohm resistance. We suggest to use this ground electrode DRY. The connection is made with a snap button and a lead wire is included in the set. Available in three different sizes.

Cod. **GDRGP02526** 2,5x25cm with press-stud connection and connection cable with Touch proof connector
Cod. **GDRGP04026** 2,5x40cm with press-stud connection and connection cable with Touch proof connector
Cod. **GDRGP07026** 2,5x70cm with press-stud connection and connection cable with Touch proof connector

Elastic Ground Electrodes

100cm long lead wire



The ground electrodes are made with a conductive elastic tissue with around 50 Ohm. It can be used dry. The connection is made with a snap button and a lead wire is included in the set. Available in three different sizes.

Cod. **GDRG002526** 2,5x25cm with press-stud connection and connection cable with Touch proof connector
Cod. **GDRG004026** 2,5x40cm with press-stud connection and connection cable with Touch proof connector
Cod. **GDRG007026** 2,5x70cm with press-stud connection and connection cable with Touch proof connector

Disposable ground electrode with


adhesive and conductive solid gel



Cod. **DENIG10026**
35x45mm disposable ground electrodes with 100cm lead wire and Touch proof connector
Cod. **DENIG15026**
35x45mm disposable ground electrodes with 150cm lead wire and Touch proof connector

Disposable ground electrode

with STUD connector



Cod. **CDES003545**
35x45mm disposable ground electrodes with stud connector

Cod. **LEAD1026300**
connection cable, 100cm, green colour with Touch proof connector

Cod. **LEDCDE0026**
new type 360° lope connection cable, 100cm with Touch proof connector