

## MYOLINE: disposable and reusable electromyography concentric Needles

Why to choos Myoline:

- A new sharper tip provides comfort for the patient
- A recorder area mark indicates the direction for insertion
- A comfortable needle holder provides a perfect handgrip
- A concentric connection allows shielding in the handgrip
- 11 sizes are available

Disposable and reusable concentric needle. Reusable EMG cable.



### ETO sterilized disposable concentric needle

Cod. D0390302511	diam. 0,30x25 PURPLE
Cod. D0390351209	diam. 0,35x12 PINK
Cod. D039035205	diam. 0,35x20 WHITE
Cod. D039035302	diam. 0,35x30 RED
Cod. D039035408	diam. 0,35x40 GREY
Cod. D039035503	diam. 0,35x50 GREEN
Cod. D0390453010	diam. 0,45x30 ORANGE
Cod. D039045401	diam. 0,45x40 BLACK
Cod. D039045504	diam. 0,45x50 YELLOW
Cod. D039050607	diam. 0,50x60 BLUE
Cod. D039060706	diam. 0,60x70 BROWN

### Non sterile Reusable concentric needle

Cod. R039035205	diam. 0,35x20 WHITE
Cod. R039035302	diam. 0,35x30 RED
Cod. R039045401	diam. 0,45x40 BLACK
Cod. R039045504	diam. 0,45x50 YELLOW
Cod. R039060706	diam. 0,60x70 BROWN

### Connection leadwire

Cod. 100000NH	Reusable EMG cable with 5-pole DIN45322 connector
Cod. 1000SCNH	Reusable EMG cable with 3 safety connector
Cod. 100200NH	Reusable EMG cable 200cm long with 5-pole DIN-45322 connector

### Other connectors available on request

## Myoline "SAFETY NEEDLE".

### The only needle which takes care of the operator

This safety needle avoids the likelihood of the operator creating a puncture as the needle and cable are molded together as one piece. The needle, cable and connector are all destroyed after use.

The only needle in the market completely sterile also in the connection leadwire.

Other sizes available on request.

### Disposable concentric needle electrode with moulded lead wire. Box 25 pcs




Cod. DC39045401

diam. 0,45x40mm ETO sterilized disposable concentric needle with 100cm long lead wire and 5-pole connector

### Disposable single fiber EMG needle

The recording surface is presented through a side port opposite to the cutting edge bevel, 5mm from the needle tip. It uses the same connector as the standard Myoline concentric needles.



Cod. DSN380453707

diam 0,45x37mm BLUE

### Disposable Teflon coated MONOPOLAR needles

An easy hand grip in a round shape is moulded on the needle and wire. A low friction Teflon coat covers the needle with the exception of recording area.

### Disposable Teflon coated monopolar needles



Cod. TFDN261500

disposable 0,35x15mm long teflon insulated monopolar needle, with 120cm long cable and Touch proof connector

Cod. TFDN263500

disposable 0,35x35mm long teflon insulated monopolar needle, with 120cm long cable and Touch proof connector


Cod. TFDN264500

disposable 0,35x45mm long teflon insulated monopolar needle, with 120cm long cable and Touch proof connector

## Disposable needles for Botulinum Toxin

This specially designed shape allows the needle to be inserted into a standard hypodermic disposable syringe or into a LUER-LOK syringe.

As the solution is injected the EMG or electric activity are recorded. The needle is coated with low friction PTFE and special care is taken to sharpen the point. Coloured lead wires are used for easy identification of the length and diameter of the needle.



Cod. **MIOBOT2530**

Disposable ETO sterilized diam. 0,30mm x 25mm needles for Botulinum Toxin blu

Cod. **MIOBOT2540**

Disposable ETO sterilized diam. 0,40mm x 25mm needles for Botulinum Toxin green

Cod. **MIOBOT3745**

Disposable ETO sterilized diam. 0,45mm x 37mm needles for Botulinum Toxin pink

Cod. **MIOBOT5050**

Disposable ETO sterilized diam. 0,50mm x 50mm needles for Botulinum Toxin white

Cod. **MIOBOT7570**

Disposable ETO sterilized diam. 0,70mm x 75mm needles for Botulinum Toxin red



0476