

DISPOCAP



A new concept in EEG recording: The DISPOCAP.

The Dispocap is a disposable cap with incorporated electrodes. The cap is knitted with a special elastic fabric in an ergonomic sheet to fit a wide range of head shape and size. The net is attached to a soft elastic with a silicon band fixed inside to guarantee the cap remains fixed in place during the recording. The electrodes are made of two parts, which are screwed together one, each side of the net. The position of the electrodes be changed if necessary. As shown the coloured wiring is placed on the outside of the cap.

The DISPOCAP is available in adult sizes: SMALL, MEDIUM, LARGE and infant MEDIUM and SMALL sizes and is designed to fit a wide range of head shape. The DISPOCAP is made up of 20 tin electrodes + 1 spare electrode to be used as REF

FOR SINGLE USE ONLY. Can be used in ICU and emergency **infective** to avoid any cross contamination.

cod DCAPSOLAC2026	Large	Adult size
cod DCAPSOMAC2026	Medium	Adult size
cod DCAPSOSAC2026	Small	Adult size
cod DCAPSOMIC2026	Medium	Infant size
cod DCAPSOSIC2026	Small	Infant size



The **DISPOCAP KIT** contains all the necessary components to record an EEG:

- 1 DISPOCAP
- 1 syringe
- 1 50gr bottle of Neurgel
- 1 spare electrodes
- 2 blunted needles

cod KDCAPSOLAC2026	Large	Adult size
cod KDCAPSOMAC2026	Medium	Adult size
cod KDCAPSOSAC2026	Small	Adult size
cod KDCAPSOMIC2026	Medium	Infant size
cod KDCAPSOSIC2026	Small	Infant size
