

Respiration sensor for sleep labs

We offer this disposable respiration sensor as an alternative to the reusable respiration belt. This sensor is embedded in an adhesive solid gel support, which guarantees an excellent, comfortable bond with the patient. This sensor is especially recommended for use during sleep recording as it does not require a chest belt and is of greater comfort to the patient.

- Peel the away the silicon paper support and attach the disposable sensor to the required part of the body
- Insert the two T.P. colour coded connectors of the disposable respiration sensor into the reusable unit
- Start the recording

DISPOSABLE Respiration sensor




Cod. RESPA000200 adult size
Cod. RESPP000200 paediatric and neonatal size



DISPOSABLE respiration sensor

CONNECTION REUSABLE UNIT



Cod. BMDS00260R
Reusable connection device for disposable respiration sensor



REUSABLE Respiration sensor



Cod. RESPA00000
Respiratory belt adult size with 2m long cables and Touch prod connector

Cod. RESPP00000
Respiratory belt paediatric size with 2m long cables and Touch prod connector



Reusable Sleep Lab Sensors

Sensors for recording body signals, respiratory airflow, limb tremors, eye movements, body position, during sleep disorder tests.

Airflow sensors are made to record the breathing signals



Cod. AIRFA00300 Reusable adult three ways respiration sensor
Cod. AIRFA00200 Reusable adult two ways respiration sensor
Cod. AIRFA00100 Reusable adult one way respiration sensor

Cod. AIRFP00300 Reusable infant three ways respiration sensor
Cod. AIRFP00200 Reusable infant two ways respiration sensor
Cod. AIRFP00100 Reusable infant one way respiration sensor

CE

Tremor movement sensor recorder



CE

Cod. TREM00000
Tremor movement sensor recorder with 2m long cables and Touch proof connector

Snoring sensor



CE

Cod. SNOR00000
This sensor is designed to record the patient's snoring vibrations he usually produces during the sleep. 2m long cable and Touch proof connector


Limb movement sensor



CE

Cod. LIMB00000
Limb movement sensor to record the periodic limb movement signals with 2m long cables and Touch proof connector.

Eye movement sensor



CE

Cod. EYE00000
This sensor is designed to record the patient's eyes movements during sleep. 2 m long cable and Touch proof connector.